

Photon-number and timing resolution of a near-IR continuous-wave source with a transition edge sensor IPS 2017

Jianwei Lee*, Lijiong Shen, Brenda Chng, Alessandro Cere,
Christian Kurtsiefer



Motivation

Measure 2nd-order correlation function with a single detector

Avoids spatial multiplexing of Hanbury-Brown-Twiss

Applicable to Satellites

Outline

Transition-Edge Sensor & Setup

Photon-Number Resolution

Time-of-Arrival Estimation


2nd-Order Correlation of Continuous Source

Transition-Edge Sensor

Photon-Number Resolving



Near unit Efficiency

Long Recovery Time ($\approx 2\mu s$) limits flux detection rates





Transition-Edge Sensor

Photon absorbed \rightarrow Temperature increases \rightarrow Resistance increases



Photon Number Discrimination


Light source: Pulsed Laser Diode
Pulse Area \propto Photon Number



Pulse Identification for Continuous Source


Discard Incomplete Pulses at Edges

Limit Background Noise




Photon Number Resolution (Continuous Source)

NO Pulse Identification




Pulse Identification



Time-of-Arrival Estimation


$$\text{Model: } A_1 S(t - t_1) + A_2 S(t - t_2)$$

1. Select 2-photon signals
(pulse area)
2. Detect Pulse edges (filter +
differentiate)
3. 2 Edges
 - Initialise using Pulse
Edges Timings
 - Least Squares Fit
4. 1 Edge
 - Initialise using Pulse
Region
 - Monte-Carlo Markov
Chain Fit




Single-Photon Detection Model

Select 5000 single-photon pulses
Correct for Vertical & Horizontal Offset
Average




Self-Correlation $G^{(2)}$

CW Laser Diode



Self-Correlation $G^{(2)}$


CW Laser Diode




Self-Correlation $G^{(2)}$ Fit Accuracy

Light Source: Laser Diode Pulse Pairs

Seperation > 200ns



Seperation = 200ns



Conclusion

Timing Resolution ≈ 250 ns

Photon Number Resolution

- $n = 0, 1$ misidentification negligible
- 1-photon event misidentified as 2-photon event: 0.8%

Applications

- Correlation measurements without spatial-mode multiplexing
- Increased Photon Flux